

Product Datasheet

Biotinylated Recombinant Human Angiopoietin-like 3 Protein, Avi His Tag (orb1550378)

Catalog Number	orb1550378
Description	Biotinylated Recombinant Human Angiopoietin-like 3 Protein, Avi His Tag
Reactivity	Human
Endotoxins	> 95% by Tris-Bis PAGE
Target	Angiopoietinlike 3
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	Biotinylated Recombinant Human Angiopoietin-like 3 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ser17-Lys219) of Human Angiopoietin-like 3 fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)