

Product Datasheet

Recombinant Mouse IL-4RA/IL-4 R alpha Protein, Fc His Tag (orb1550765)

Catalog Number	orb1550765
Category	Proteins
Description	Recombinant Mouse IL-4RA/IL-4 R alpha Protein, Fc His Tag
Target	IL4RA
Reactivity	Mouse
Purity	> 90% by SDS-PAGE.
Expression System	HEK293 Cells
Endotoxins	> 90% by SDS-PAGE.
Source	Recombinant Mouse IL-4RA/IL-4 R alpha Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ile26-Arg233) of mouse IL-4RA/IL-4 R alpha (Accession #P16382) fused with an Fc tag at the C-terminus.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)