

Product Datasheet

Recombinant Mouse IL-4 Protein, His Tag (orb1551183)

Catalog Number	orb1551183
Description	Recombinant Mouse IL-4 Protein, His Tag
Reactivity	Mouse
Endotoxins	> 95% by SDS-PAGE.
Target	IL4
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	Recombinant Mouse IL-4 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (His Tag23-Ser140) of mouse IL-4 (Accession #NP_067258.1) fused with a 6His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)