

## Product Datasheet

### Recombinant Human ISG15 Protein (orb1551428)

<b>Catalog Number</b>	orb1551428
<b>Category</b>	Proteins
<b>Description</b>	Recombinant Human ISG15 Protein
<b>Target</b>	ISG15
<b>Reactivity</b>	Human
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Expression System</b>	E. coli
<b>Endotoxins</b>	> 97% by SDS-PAGE.
<b>Source</b>	Recombinant Human ISG15 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gly2-Gly157) of human ISG15 (Accession #NP_005092.1).
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)