

Product Datasheet

Recombinant Human AHSG Protein, His Tag (orb1553085)

Catalog Number	orb1553085
Category	Proteins
Description	Recombinant Human AHSG Protein, His Tag
Target	AHSG
Reactivity	Human
Purity	≥90% as determined by SDS-PAGE
UniProt ID	P02765
Expression System	E. coli
Endotoxins	≥90% as determined by SDS-PAGE
Source	This protein contains the human AHSG(Ala19-Val367) was fused with the N-terminal His Tag and expressed in E. coli.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)