

Product Datasheet

MYH10 Rabbit Polyclonal Antibody (orb1562744)

Catalog Number	orb1562744
Category	Antibodies
Description	Purified Rabbit Anti Human MYH10 Polyclonal antibody.
Target	MYH10
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ . Or dissolved in sterile H ₂ O with 50% glycerol.
Purification	Purified by antigen affinity column.
Immunogen	Recombinant fragment corresponding to 1507-1606aa of human MYH10.
MW	230KD
Tested applications	ELISA, WB
Dilution range	WB_value: 1:800-1:2,000. ELISA_value: 1:20,000-1:80,000
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Biorbyt Ltd.	Biorbyt LLC
7 Signet Court, Swann Road	68 TW Alexander Drive
Cambridge	Research Triangle Park
CB5 8LA	Durham
United Kingdom	NC 27713
Email: info@biorbyt.com , support@biorbyt.com	United States
Phone: +44 (0)1223 859353 Fax: +1 (415) 651-8558	Email: info@biorbyt.com , support@biorbyt.com
	Phone: +1 (415) 906-5211 Fax: +1 (415) 651-8558

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)