

Product Datasheet

HA Tag (3C6) Mouse mAb Antibody (orb1562813)

Catalog Number	orb1562813
Category	Antibodies
Description	HA Tag (3C6) Mouse mAb Antibody
Clonality	Monoclonal
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Other
Form/Appearance	Liquid
Concentration	1.48 mg/ml
Buffer/Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified
Immunogen	HA synthetic peptide[YPYDVPDYA]conjugated to KLH.
Tested applications	IP, WB
Dilution range	WB: 1/1000-1/10000 IP: 1/20
Clone Number	B5E8XH0X10
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.



Western blot analysis of HA Tag (3C6) in CHO-K1 (C, F) or CHO-K1 transfected with different HA fusion proteins (A, B, D, E), using HATag antibody.



Immunoprecipitation analysis of HA Tag (3C6) in HeLa cell lysates transfected with HA tagged fusion protein using HATag antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558