

Product Datasheet

GFP Tag Mouse mAb (orb1562835)

Catalog Number	orb1562835
Category	Antibodies
Description	GFP Tag Mouse mAb
Target	GFP
Clonality	Monoclonal
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Other
Form/Appearance	Liquid
Concentration	1 mg/ml
Buffer/Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified
Immunogen	Recombinant Protein of GFP
UniProt ID	Q963I9
MW	Calculated MW: 27 kDa; Observed MW:
Tested applications	IP, WB
Dilution range	WB: 1/1000-1/10000 IP: 1/20

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558


Clone Number 6A2-1H8-5G7

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only



1) Input (control) 2) IP products, antibody dilution 1:200.



Western blot analysis of GFPtag in CHO-K1 transfected by different GFPfusion proteins (A, B, C) or nontransfected (D) cell lysates using GFPtag antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558