

Product Datasheet

GFP/eGFP Tag Mouse mAb (orb1562842)


Description	GFP/eGFP Tag Mouse mAb
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant EYFP.tag full length expressed in E.coli.
Target	eGFP
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	2.28 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	3F11-B12-E7-F12
MW	Calculated MW: 27 kDa; Observed MW:
Uniprot ID	C8CHS1
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Western blot analysis of EGFP/EYFP in Rosetta E.coli induced by 0 (A) or 0.1mM (B) IPTG using EGFP/EYFP antibody.



Western blot analysis of EGFP/EYFP in CHO-K1 transfected by different EGFPfusion proteins or nontransfected (A) cell lysates using EGFP/EYFP antibody.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)