

Product Datasheet

Human FGL2 Protein, His Tag (orb1570937)

Catalog Number	orb1570937
Category	Proteins
Description	Human FGL2 (204-439) Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Pro 204 - Pro 439 (Accession # Q14314-1).
Target	FGL2
Conjugation	Unconjugated
Tag	N-10×His
Reactivity	Human
Form/Appearance	Powder
Buffer/Preservatives	PBS, pH 7.4
Purity	90%
UniProt ID	Q14314
MW	29.6 kDa
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 29.6 kDa. The protein migrates as 36-41 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Expression System	HEK293
Endotoxins	1.0 EU per µg
Expression Region	Pro 204 - Pro 439

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)