

Product Datasheet

Phospho-GSK-3 Beta (Ser9) Rabbit Polyclonal Antibody (PerCP-Cy7) (orb1579646)

Description	Phospho-GSK-3 Beta (Ser9) Rabbit Polyclonal Antibody (PerCP-Cy7)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PerCP/Cy7
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Ser9 TT(p-S)FA
Target	GSK3B
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin and 0.02% Proclin300.
Concentration	1mg/ml
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	Flow-Cyt=1ug/Test, ICC/IF=1:100, IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)