

Product Datasheet

Phospho-RAD51 (Thr309) Rabbit Polyclonal Antibody (PerCP-Cy5.5) (orb1582742)

Catalog Number	orb1582742
Category	Antibodies
Description	Phospho-RAD51 (Thr309) Rabbit Polyclonal Antibody (PerCP-Cy5.5) is a PerCP/Cy5.5 conjugated antibody targeting RAD51. This antibody is suitable for FC, ICC, IF. It exhibits reactivity with Human, Mouse samples.
Target	RAD51
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PerCP/Cy5.5
Reactivity	Human, Mouse
Predicted Reactivity	Bovine, Canine, Gallus, Porcine, Rabbit, Rat, Sheep
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin and 0.02% Proclin300.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RAD51 around the phosphorylation site of Thr309 GE(p-T)RI
Tested applications	FC, ICC, IF

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	Flow-Cyt=1ug/Test, ICC/IF=1:100-500, IF=1:100-500
Antibody Type	Primary Antibody
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)