

Product Datasheet

Alpha smooth muscle Actin Recombinant Rabbit Monoclonal Antibody (BF680) (orb1586527)

Catalog Number	orb1586527
Category	Antibodies
Description	Alpha smooth muscle Actin Recombinant Rabbit Monoclonal Antibody (BF680) is a BF680 conjugated antibody targeting ACTA2. This antibody is suitable for ICC, IF. It exhibits reactivity with Human, Mouse, Rat samples.
Target	ACTA2
Clonality	Recombinant
Species/Host	Rabbit
Isotype	IgG
Conjugation	BF680
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Zebrafish
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthetic peptide derived from human Actin alpha
Tested applications	FC, ICC, IF
Dilution range	ICC/IF=1:50-200, IF=1:100-500, Flow-Cyt=1ug/Test

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Antibody Type	Primary Antibody
Clone Number	B4D32
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)