

Product Datasheet

Phospho-IGF1R (Tyr980) Rabbit Polyclonal Antibody (PerCP-Cy7) (orb1591658)

| | |
|-----------------------------|---|
| Catalog Number | orb1591658 |
| Category | Antibodies |
| Description | Phospho-IGF1R (Tyr980) Rabbit Polyclonal Antibody (PerCP-Cy7) is a PerCP/Cy7 conjugated antibody targeting IGF1R. This antibody is suitable for FC, IF. It exhibits reactivity with Human, Mouse samples. |
| Target | IGF1R |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | IgG |
| Conjugation | PerCP/Cy7 |
| Reactivity | Human, Mouse |
| Predicted Reactivity | Rat |
| Concentration | 1mg/ml |
| Buffer/Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin and 0.02% Proclin300. |
| Purification | Affinity purified by Protein A |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human IGF1R around the phosphorylation site of Tyr980 PE(p-Y)FS |
| Tested applications | FC, IF |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Dilution range | Flow-Cyt=1µg/Test, IF=1:100-500 |
| Antibody Type | Primary Antibody |
| Storage | Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)