

Product Datasheet

Phospho-NFkB1 (Ser373) Rabbit Polyclonal Antibody (APC-Cy7) (orb1592772)

Catalog Number	orb1592772
Category	Antibodies
Description	Phospho-NFkB1 (Ser373) Rabbit Polyclonal Antibody (APC-Cy7)
Target	NFkB1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	APC/Cy7
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Equine, Porcine, Sheep
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NF KappaB p105 around the phosphorylation site of Ser373 SD(p-S)FG
Tested applications	FC, IF
Dilution range	Flow-Cyt=1µg/Test, IF=1:100-500

Antibody Type	Primary Antibody
----------------------	------------------

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)