

Product Datasheet

mScarlet Mouse Monoclonal Antibody (orb1597254)

Catalog Number	orb1597254
Category	Antibodies
Description	mScarlet Mouse Monoclonal Antibody
Target	mScarlet
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Other
Predicted Reactivity	Other
Form/Appearance	Liquid
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	Recombinant mScarlet proteins
MW	26.4 kDa
Tested applications	ELISA, FC, ICC, IF, WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution range	WB: 1:2000-5000; IF/ICC: 1:100-1:500; FC: 2ug/test
Antibody Type	Primary Antibody
Clone Number	B32F9
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T-mScarlet cells using mScarlet Mouse Monoclonal Antibody orb1597254 with dilution at 1:3000.



Immunofluorescent analysis of 293T-mScarlet cells using mScarlet Mouse Monoclonal Antibody orb1597254 with dilution at 1:100.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Immunofluorescent analysis of 293T-mScarlet cells using mScarlet Mouse Monoclonal Antibody orb1597254 with dilution at 1:100.



Immunofluorescent analysis of 293T-mScarlet cells using mScarlet Mouse Monoclonal Antibody orb1597254 with dilution at 1:100.



Flow Cytometry analysis of 293t-mScarlet cells (1×10^5 cells) using orb1597254 mScarlet Mouse Monoclonal Antibody with 2ug /test (green), blank (grey), isotype control (blue).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558