

Product Datasheet

Phospho-IKB alpha (Ser32 + Ser36) Rabbit Polyclonal Antibody (PerCP-Cy5.5) (orb1603004)

Catalog Number	orb1603004
Category	Antibodies
Description	Phospho-IKB alpha (Ser32 + Ser36) Rabbit Polyclonal Antibody (PerCP-Cy5.5)
Target	NFKBIA
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PerCP/Cy5.5
Reactivity	Human, Rat
Predicted Reactivity	Bovine, Canine, Gallus, Mouse, Porcine, Rabbit, Sheep
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin and 0.02% Proclin300.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IKB alpha around the phosphorylation site of Ser32/36 HD(p-S)GLD(p-S)M
Tested applications	ICC, IF
Dilution range	ICC/IF=1:100-500, IF=1:100-500
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Storage Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)