

Product Datasheet

Phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199) Rabbit Polyclonal Antibody (APC-Cy5.5) (orb1605197)

Description	Phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199) Rabbit Polyclonal Antibody (APC-Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	APC/Cy5.5
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467 RL(p-Y)EE
Target	PIK3R1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	Flow-Cyt=1ug/test, IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)