


Product Datasheet

Rat IgG (orb1610937)

Catalog Number	orb1610937
Category	Proteins
Description	This product is Rat IgG, also known as Rat Immunoglobulin G, native; Native Rat IgG. It was purified by affinity purified by Protein A and provided in lyophilized form.
Form/Appearance	Lyophilized
Buffer/Preservatives	0.01M PBS (pH7.4).
Purification	affinity purified by Protein A
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558