

Product Datasheet

Mouse IgG (orb1611040)

Catalog Number	orb1611040
Category	Proteins
Description	This product is Mouse IgG, also known as Native Mouse IgG. It is provided as a lyophilized form (buffer before lyophilization is 0.01M PBS, pH 7.4), affinity purified by Protein G. Store at -20°C. Avoid repeated freeze/thaw cycles.
Form/Appearance	Lyophilized
Buffer/Preservatives	0.01M PBS (pH7.4).
Purification	affinity purified by Protein G
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)