

## Product Datasheet

### Human IgM (orb1611042)

|                             |                                                                 |
|-----------------------------|-----------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1611042                                                      |
| <b>Category</b>             | Proteins                                                        |
| <b>Description</b>          | Human IgM Antibody                                              |
| <b>Form/Appearance</b>      | Lyophilized or Liquid                                           |
| <b>Buffer/Preservatives</b> | 0.01M PBS (pH7.4).                                              |
| <b>MW</b>                   | 950 kDa                                                         |
| <b>Storage</b>              | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                 | For research use only                                           |
| <b>Expiration Date</b>      | 12 months from date of receipt.                                 |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)