

Product Datasheet

RSK1 p90(Phospho-T359+S363) Conjugated Antibody (orb1612518)

Description	RSK1 p90(Phospho-T359+S363) Conjugated Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Antibody Type	Primary Antibody
Source	Recombinant Rabbit Monoclonal Ab
Uniprot ID	Q9BXW9
Dilution Range	working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)