

Product Datasheet

JNK1/2/3 (Phospho-Tyr185) Polyclonal Polyclonal Conjugated Antibody (orb1630364)

Catalog Number	orb1630364
Category	Antibodies
Description	JNK1/2/3 (Phospho-Tyr185) Polyclonal Polyclonal Conjugated Antibody
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Reactivity	Human
UniProt ID	P53779, P45983, P45984
Dilution range	working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Antibody Type	Primary Antibody
Source	Rabbit Polyclonal Ab
Storage	Store at 4°C for 6 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)