

Product Datasheet

Amyloid β A4 (Phospho-Thr743/668) Polyclonal Conjugated Antibody (orb1630875)

Catalog Number	orb1630875
Category	Antibodies
Description	Amyloid β A4 (Phospho-Thr743/668) Polyclonal Conjugated Antibody
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Reactivity	Human, Mouse, Rat
UniProt ID	P05067
Dilution range	working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Antibody Type	Primary Antibody
Source	Rabbit Polyclonal Ab
Storage	Store at 4°C for 6 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)