

## Product Datasheet

# MKK3 (Phospho-Ser189) Polyclonal Conjugated Antibody (orb1631044)

<b>Catalog Number</b>	orb1631044
<b>Category</b>	Antibodies
<b>Description</b>	MKK3 (Phospho-Ser189) Polyclonal Conjugated Antibody
<b>Clonality</b>	Polyclonal
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Reactivity</b>	Human, Mouse, Rat
<b>UniProt ID</b>	<b>P46734</b>
<b>Dilution range</b>	working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
<b>Antibody Type</b>	Primary Antibody
<b>Source</b>	Rabbit Polyclonal Ab
<b>Storage</b>	Store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)