

Product Datasheet

G3BP-1 (Phospho-Ser232) Polyclonal Conjugated Antibody (orb1631087)

Catalog Number	orb1631087
Description	G3BP-1 (Phospho-Ser232) Polyclonal Conjugated Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Polyclonal
Antibody Type	Primary Antibody
Source	Rabbit Polyclonal Ab
Uniprot ID	Q13283
Dilution Range	working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)