

Product Datasheet

Recombinant Mouse C-C motif chemokine 5 protein(Ccl5) (Active) (orb1650791)

Catalog Number	orb1650791
Category	Proteins
Description	Recombinant Mouse C-C motif chemokine 5 protein(Ccl5) (Active)
Tag	NO-tagged
Form/Appearance	0.2 ?m filtered PBS, pH 7.4 ,lyophilized
Purity	>97% as determined by SDS-PAGE and HPLC.
Protein Sequence	SPYGSDDTTPC CFAYLSLALP RAHVKEYFYT SSKCSNLAVV FVTRRRNRQVC ANPEKKWVQE YINYLEMS
UniProt ID	P30882
MW	7.9 kDa
Expression System	E.Coli
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human T lymphocytes is in a concentration range of 1.0-10 ng/ml.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)