

Product Datasheet

Recombinant Mouse C-C motif chemokine 24 protein(Ccl24) (Active) (orb1650807)

Catalog Number	orb1650807
Category	Proteins
Description	Recombinant Mouse C-C motif chemokine 24 protein(Ccl24) (Active)
Tag	NO-tagged
Form/Appearance	0.2 ?m filtered PBS, pH 7.4 ,lyophilized
Purity	>97% as determined by SDS-PAGE and HPLC.
Protein Sequence	VTIPSSCCTS FISKKIPENR VVSYQLANGS ICPKAGVIFI TKKGHKICTD PKLLWVQRHI QKLDAKKNQP SKGAKAVRTK FAVQRRRGNS TEV
UniProt ID	Q9JKC0
MW	10.3 kDa
Expression System	E.Coli
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using murine lymphocytes is in a concentration of 10-100 ng,ml.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)