

## Product Datasheet

# Human Phosphatidylinositol-binding clathrin assembly protein (PICALM) ELISA Kit (orb1654976)

<b>Catalog Number</b>	orb1654976
<b>Category</b>	Assays and Kits
<b>Description</b>	Human Phosphatidylinositol-binding clathrin assembly protein (PICALM) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.118 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)