

Product Datasheet

Mouse A Disintegrin And Metalloprotease 8 (ADAM8) ELISA Kit (orb1655780)

Catalog Number	orb1655780
Category	Assays and Kits
Description	Mouse A Disintegrin And Metalloprotease 8 (ADAM8) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 7.81-500 pg/mL. Sensitivity: 3.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)