

## Product Datasheet

# Human TNF receptor-associated factor 3 (TRAF3) ELISA Kit (orb1657455)

|                            |                                                                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb1657455                                                                                                                                                                                    |
| <b>Category</b>            | Assays and Kits                                                                                                                                                                               |
| <b>Description</b>         | Human TNF receptor-associated factor 3 (TRAF3) ELISA Kit                                                                                                                                      |
| <b>Reactivity</b>          | Human                                                                                                                                                                                         |
| <b>Tested applications</b> | ELISA                                                                                                                                                                                         |
| <b>Application notes</b>   | Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.3 pg/mL. Assay Time: 1-4.5 h.<br>Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum,<br>Plasma, Other biological fluids. |
| <b>Note</b>                | For research use only                                                                                                                                                                         |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                |

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)