

Product Datasheet

Mouse GREM2 ELISA Kit (orb1657745)

Catalog Number	orb1657745
Category	Assays and Kits
Description	Mouse GREM2 ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Grem2. Gene_full_name: Gremlin-2. Protein_name: Gremlin-2. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Research_category: signal transduction cytoskeleton ecm extracellular matrix structures bone stem cells embryonic stem cells cytokines. Tissue_specificity: Highly expressed in the ovary followed by brain spleen colon kidney and uterus. In ovary expressed in granulosa cells of selective early antral follicles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)