

Product Datasheet

Mouse R-Spondin-3 ELISA Kit (orb1657749)

Catalog Number	orb1657749
Category	Assays and Kits
Description	Mouse R-Spondin-3 ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Rspo3. Gene_full_name: R-spondin-3. Protein_name: R-spondin-3. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Research_category: stem cells signaling pathways wnt secreted developmental biology reproduction placental development. Tissue_specificity: Highly expressed in endothelial cells(PubMed:26766444).
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)