

## Product Datasheet

### Human Livin ELISA Kit (orb1657854)

<b>Catalog Number</b>	orb1657854
<b>Category</b>	Assays and Kits
<b>Description</b>	Human Livin ELISA Kit
<b>Reactivity</b>	Human
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detect Range: 0.78 ng/mL-50 ng/mL. Sensitivity: 20 pg/mL. Gene_name: BIRC7. Gene_full_name: Interleukin-17C. Protein_name: Baculoviral IAP repeat-containing protein 7. Sample Type: cell culture supernates serum cell lysates and tissue homogenates. Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: cell biology apoptosis intracellular survivin iaps cancer invasion microenvironment death receptors & ligands cell death receptors. Sequence_similarities: Belongs to the IAP family. Tissue_specificity: Isoform 1 and isoform 2 are expressed at very low levels or not detectable in most adult tissues. Detected in adult heart placenta lung lymph node spleen and ovary and in several carcinoma cell lines. Isoform 2 is detected in fetal kidney heart and spleen and at lower levels in adult brain skeletal muscle and peripheral blood leukocytes.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)