

## Product Datasheet

### Human CRLF2/TSLP R ELISA Kit (orb1657893)

<b>Catalog Number</b>	orb1657893
<b>Category</b>	Assays and Kits
<b>Description</b>	Human CRLF2/TSLP R ELISA Kit
<b>Reactivity</b>	Human
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CRLF2. Gene_full_name: Cytokine receptor-like factor 2. Protein_name: Cytokine receptor-like factor 2. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines. Sequence_similarities: Belongs to the type I cytokine receptor family. Type 5 subfamily. Tissue_specificity: Expressed in heart skeletal muscle kidney and adult and fetal liver. Primarily expressed in dendrites and monocytes. Weakly expressed in T-cells.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)