

Product Datasheet

Human IL-18R1 ELISA Kit (orb1657910)

Catalog Number	orb1657910
Category	Assays and Kits
Description	Human IL-18R1 ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: IL18R1. Gene_full_name: Interleukin-18 receptor 1. Protein_name: Interleukin-18 receptor 1. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. cytokines interleukins. Sequence_similarities: Belongs to the interleukin-1 receptor family. Tissue_specificity: Expressed in lung leukocytes spleen liver thymus prostate small intestine colon placenta and heart and is absent from brain skeletal muscle pancreas and kidney. High level of expression in Hodgkin disease cell lines.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)