

Product Datasheet

Mouse TNFRSF4/OX40 ELISA Kit (orb1658001)

Catalog Number	orb1658001
Category	Assays and Kits
Description	Mouse TNFRSF4/OX40 ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFRSF4. Gene_full_name: Tumor necrosis factor receptor superfamily member 4. Protein_name: Tumor necrosis factor receptor superfamily member 4. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology adaptive immunity b cells cd t cells innate immunity cytokines tnf superfamily immune system diseases autoimmune stem cells hematopoietic progenitors lymphoid t lymphocytic lineage regulatory t cells.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)