

Product Datasheet

Human CSF1R/M-CSFR ELISA Kit (orb1658114)

Catalog Number	orb1658114
Category	Assays and Kits
Description	Human CSF1R/M-CSFR ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CSF1R. Gene_full_name: Macrophage colony-stimulating factor 1 receptor. Protein_name: Macrophage colony-stimulating factor 1 receptor. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular angiogenesis cytokines immunology innate immunity macrophage inflamm. csfs signal transduction protein phosphorylation ser thr kinases kinases. Sequence_similarities: Belongs to the protein kinase superfamily. Tyr protein kinase family. CSF-1 PDGF receptor subfamily. Tissue_specificity: Expressed in bone marrow and in differentiatedblood mononuclear cells.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)