

Product Datasheet

Human IL-1RA/IL1RN ELISA Kit (orb1658129)

Catalog Number	orb1658129
Category	Assays and Kits
Description	Human IL-1RA/IL1RN ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 2 pg/mL. Gene_name: IL1RN. Gene_full_name: Interleukin-1 receptor antagonist protein. Protein_name: Interleukin-1 receptor antagonist protein. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: IL1RN. Tissue_specificity: The intracellular form of IL1RN is predominantly expressed in epithelial cells.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)