

## Product Datasheet

### Human CXCL4/PF4 ELISA Kit (orb1658166)

<b>Catalog Number</b>	orb1658166
<b>Category</b>	Assays and Kits
<b>Description</b>	Human CXCL4/PF4 ELISA Kit
<b>Reactivity</b>	Human
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 15 pg/mL. Gene_name: PF4. Gene_full_name: Platelet factor 4. Protein_name: Platelet factor 4. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Research_category: cardiovascular angiogenesis inhibitors immunology innate immunity macrophage inflamm. blood platelets chemokines alpha chemokines (cxc) cancer invasion microenvironment angiogenic inhibitory factors atherosclerosis vascular inflammation leukocyte recruitment kits lysates other kits elisa kits atherosclerotic proteins elisa kits.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)