

## Product Datasheet

### Rat MIP-3 Alpha/CCL20 ELISA Kit (orb1658275)

<b>Catalog Number</b>	orb1658275
<b>Category</b>	Assays and Kits
<b>Description</b>	Rat MIP-3 Alpha/CCL20 ELISA Kit
<b>Reactivity</b>	Rat
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detect Range: 7.8 pg/mL-500 pg/mL. Sensitivity: 1 pg/mL. Gene_name: CCL20. Gene_full_name: C-C motif chemokine 20. Protein_name: C-C motif chemokine 20. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) signal transduction signaling pathway g protein signaling gpcr kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits. Tissue_specificity: Low levels in thymus and lung.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)