

## Product Datasheet

# Rat Aifm1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit (orb1661990)

|                            |                                                                                                                                                                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb1661990                                                                                                                                                                                                                                               |
| <b>Category</b>            | Assays and Kits                                                                                                                                                                                                                                          |
| <b>Description</b>         | Rat Aifm1/Apoptosis-inducing factor 1, mitochondrial ELISA Kit                                                                                                                                                                                           |
| <b>Reactivity</b>          | Rat                                                                                                                                                                                                                                                      |
| <b>Tested applications</b> | ELISA                                                                                                                                                                                                                                                    |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL. Sensitivity: 14.3 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Aifm1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                    |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                                                           |

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)