

Product Datasheet

Mouse Adam10/Disintegrin and metalloproteinase domain-containing protein 10 ELISA Kit (orb1662065)

Catalog Number	orb1662065
Category	Assays and Kits
Description	Mouse Adam10/Disintegrin and metalloproteinase domain-containing protein 10 ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 0.78-50 ng/mL. Sensitivity: 0.27 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Adam10. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)