

Product Datasheet

Bovine CCL3/C-C motif chemokine 3 ELISA Kit (orb1662246)

Catalog Number	orb1662246
Category	Assays and Kits
Description	Bovine CCL3/C-C motif chemokine 3 ELISA Kit
Reactivity	Bovine
Tested applications	ELISA
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.8 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: CCL3. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)