

Product Datasheet

Human Ubiquitin Specific Peptidase 14 (USP14) ELISA Kit (orb1662911)

Catalog Number	orb1662911
Category	Assays and Kits
Description	Human Ubiquitin Specific Peptidase 14 (USP14) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.108 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)