

Product Datasheet

Human Proislet Amyloid Polypeptide (ProIAPP) ELISA Kit (orb1663931)

Catalog Number	orb1663931
Category	Assays and Kits
Description	Human Proislet Amyloid Polypeptide (ProIAPP) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 10.15 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)