

Product Datasheet

Human Protein Phosphatase 1, Regulatory Subunit 1B (PPP1R1B) ELISA Kit (orb1663984)

Catalog Number	orb1663984
Category	Assays and Kits
Description	Human Protein Phosphatase 1, Regulatory Subunit 1B (PPP1R1B) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.072 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)