

Product Datasheet

Human Myosin Heavy Chain 11, Smooth Muscle (MYH11) ELISA Kit (orb1664646)

Catalog Number	orb1664646
Category	Assays and Kits
Description	Human Myosin Heavy Chain 11, Smooth Muscle (MYH11) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 0.781-50 ng/mL. Sensitivity: 0.27 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)