

## Product Datasheet

# Rat Mixed Lineage Kinase Domain Like Protein (MLKL) ELISA Kit (orb1664814)

|                            |   |
|----------------------------|---|
| <b>Catalog Number</b>      | orb1664814  |
| <b>Category</b>            | Assays and Kits   |
| <b>Description</b>         | Rat Mixed Lineage Kinase Domain Like Protein (MLKL) ELISA Kit   |
| <b>Reactivity</b>          | Rat   |
| <b>Tested applications</b> | ELISA   |
| <b>Application notes</b>   | Detect Range: 31.2-2000 pg/mL. Sensitivity: 13.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids. |
| <b>Note</b>                | For research use only   |
| <b>Expiration Date</b>     | 6 months from date of receipt.  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)