

Product Datasheet

General Leukotriene B4 (LTB4) ELISA Kit (orb1664979)

Catalog Number	orb1664979
Description	General Leukotriene B4 (LTB4) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.28 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)